

Argon in electrosurgery

argon in electrosurgery

efficiency oriented performance

Argon plasma coagulation ensures fast and efficient coagulation of large, heavily bleeding surfaces. Provides effectual devitalisation of tissues, e.g. neoplastic tumours. Argon plasma coagulation means less blood loss and less tissue damage. With penetration depth limited to 3 mm, it is particularly recommended for areas of high perforation risk.



argon in electrosurgery

confidence through technology





Argon coagulation uses the phenomenon of good condution of high frequency current by ionised argon. Argon is a chemically inert gas, devoid of physiological effects and non combustible. Under the effect of current, it becomes ionised and forms a plasma cloud in which electric arcs are formed.

In argon coagulation, there is no contact of the active electrode with the tissue, and the distance between the surgical instrument and the tissue in open surgery is up to about 5 mm, and in endoscopic surgery up to about 3 mm. The thermal effect occurs at the time when a spark jumps from the active electrode tip to the tissue. The length of the plasma arc between the probe tip and the tissue depends on the selected power, resistance of the target tissue and argon flow rate. Usually the distance between the active electrode and the tissue is 3 to 5 mm, depending on the selected coagulation parameters.

When using argon coagulation, observe all precautions for standard monopolar coagulation. Read the instructions for using argon components. Class 4.8 (99.998%) or 5.0 (99.999%) argon is used for argon coagulation.



argon safety features

confidence through safety

The main advantage of argon coagulation is constant, minimum depth of the thermal effect. Owing to the limited depth of tissue damage during argon plasma coagulation, the risk of perforation is minimised; therefore, this method can be safely used in thin walled organs.

In the case of classical contact electrocoagulation, the thermal effect reaches deaper into the tissue; this is associated with a risk of gastrointestinal perforation. In argon coagulation, the plasma arc occurs in the tissues that have the lowest electrical resistance. The tissue through which the current has flowed achieves rapid haemostasis and as a result its electrical resistance increases. It means that at that site electric arcs will not form any more, so the coagulation depth will not increase and it will be maintained within the limits of 2 to 3 mm.

As argon plasma is a good conductor, the desired effect is obtained with significantly less power compared to standard high voltage coagulation and the amount of heat delivered to the patient's tissues is lower.

advantages of argon

setting performance goals

Efficiency enhancing features of argon plasma coagulation:

- immediate hemostasis helps efficiently coagulate large areas of bleeding surface
- penetration depth limited to approximately 3 mm minimizes risk of perforation
- tissue carbonization is minimal compared to standard electrocoagulation
- no tissue vaporization minimizes the risk of perforation
- no contact between the applicator and tissue means no tissue adhesion
- · less surgical smoke gives good visibility of operating area
- · reduced smoke eliminates unpleasant odors
- · precise application of thermal energy results in reducing procedure time



ARGON COAG

The mode is used for non-contact coagulation of the surface of bleeding tissues. It eliminates smoke and smell. It ensures very shallow and gentle coagulation. Instruments: rigid argon electrodes for coagulation.



ENDO ARGON

Argon-enhanced monopolar coagulation for endoscopic procedures. It ensures very shallow and gentle coagulation. It is necessary when there is a risk of perforation. The absence of smoke ensures the perfect visibility of the operative field. Instruments: flexible argon probes.



PULSE ARGON

Argon-enhanced pulsed monopolar coagulation. It is used in gastroenterology to control bleeding. It enables the precise delivery of energy doses exactly to the bleeding site. Instruments: flexible argon probes.



ARGON CUT

Argon-enhanced monopolar cutting. The use of argon reduces the amount of smoke and smell. The thermal damage to the tissues is reduced and bleeding control is improved. This function is particularly desirable during procedures that require intensive use of the unit.

Instruments: needle or lancet type argon electrodes.

applications

Argon coagulation has broad application in local treatment of cancer, both in the case of advanced tumour resection and in the treatment of benign or pre-cancerous lesions.

The method of coagulation in argon plasma is equally effective when treating small adenomas of the large intestine, and, first of all, in tumour resection at locations where the risk of perforation is high.

Due to its convenience and safety, argon coagulation is often used during procedures performed for non oncological indications - endoscopic bleeding control and destruction of vascular lesions in the gastrointestinal tract.

General Surgery	open liver surgery, e.g. superficial bleeding after partial hepatectomy	
	 abdominal surgery 	
	 breast surgery, e.g. breast reconstruction, breast reduction, removing breast tumours 	
Bronchoscopy	 superficial bleeding 	
	 benign endobronchial tumors as papillomatoses, granuloma, lipoma, hemangiomas 	
	 recanalization of malign stenoses of the respiratory tract 	
	 stent ingrowth/overgrowth 	
	• scar stenoses	
	post-interventional conditioning of resected area, e.g. after cryosurgery	
Pulmonology	hemoptysis	
	surface haemorrhages	
	 benign endobronchial tumors (eg. papillomatosis, granulomas polyps in the trachea, lipomas, hemangiomas) 	
	 recanalization of malignant stenoses of the respiratory tract 	
	 stent ingrowth / overgrowth 	
	scar stenoses	
Gastroenterology	 haemorrhage from angiodysplastic lesions 	
	 haemorrage from polypectomy sites 	
	 devitalization of remaining tissue after polypectomy 	
	 erosions or ulcers or oozing of blood due to vascular penetration by tumors 	
	 residual sessile adenoma tissue 	
	 stenosing tumors 	
	 stent ingrowth 	
	• colitis	
	 bleeding gastric or colon carcinoma 	
	watermelon stomach	
Otolaryngology	tonsillectomy	
	 therapy of subglottic and tracheal lesions (e.g. granulomas after laser surgery, papillomatosis, subglottic stenosis) 	
	 microsurgery of soft palate (e.g. sleep apnea, fibromas) 	
	 therapy of superficial muscoal lesions (e.g. leucoplakia, hemangiomas, granulomas, precancerosis) 	
	 applications in the nasal cavity (e.g. nasal hemorrhaging, hyperplasia of the nasal concha) 	
Gynecology	treating the uterus during a myomectomy	
-,01	laparoscopic surgery for endometriosis	

instruments for endoscopic and open surgical procedures

setting the stage

Our offer includes a comprehensive range of electrosurgical products for argon plasma coagulation, completely equipped with accessories and surgical instruments - for both endoscopic procedures and open surgery. We also offer an argon tip for laparoscopic applications.

All instruments for argon coagulation offered by EMED are intended for multiple use. They can be sterilised in an autoclave at 134°C.

TROLLEYS FOR ARGON PLASMA GENERATORS

080-100

SpectrumLine trolley with argon cylinder case for electrosurgical units







electrosurgical units

with argon module



electrosurgical instruments

SDS argon instruments for open surgery and laparoscopy



electrosurgical instruments

SDS argon endoscopic instruments

	432-46A	Monopolar cable for argon flexible electrode, flat connector, L: 3.5m, SDSA plug
113-	432-46S	Monopolar cable for argon flexible electrode, L: 3.5m, SDS/LuerLock plug, flat connector
	932-148	Flexible argon probe, reusable, TBS, dia. 1.5mm, length 1.5m
	932-149	Flexible argon probe, reusable, GIT, dia. 2.3mm, lenght 2.2m
	932-150	Flexible argon probe, reusable, GIT, dia. 3.2mm, length 2.2m
	932-151	Flexible argon probe, reusable, TBS, dia. 1.5mm, length 3m
flat type plug / złącze płaskie	932-152	Flexible argon probe, reusable, GIT, dia. 2.3mm, length 3m

electrosurgical instruments

argon endoscopic instruments



argon accessories



contact us





EMED products are available all over the world. See **www.emed.pl** for contact details.

EMED SP. Z O. O. SP. K.

Ryżowa 69a, 05-816 Opacz-Kolonia Poland tel: + 48 22 723 08 00 export@emed.pl www.emed.pl





Attention! This brochure does not replace instructions for use! Refer to instructions for use!

Copyright© EMED. 2274 All rights reserved. Any copying, distribution, publishing in whole or in part without written EMED permission is prohibited.

REF NO 2021.02.K03.EN