

# Electrosurgery

# confidence in performance

As advanced new technologies are developed, so there is constant evolution in surgical techniques, leading to dramatic advances in every area of surgical procedure. In particular these changes apply to electrosurgery.

# confidence in performance

## Tradition and experience

EMED has been present in the market since 1995. From the beginning, we have operated in the medical segment. Electrosurgery is our domain. Starting from a single model of electrosurgical unit which entered into production in 1997, we gradually extended our offer. Today, we offer a whole range of electrosurgical units fully equipped with accessories and instruments. Our electrosurgical products vary from small units for use in outpatient settings to comprehensive electrosurgical systems.

## Our mission

Based on our long standing experience in electrosurgical units' production all together with our customer's knowledge, we want to ensure constant development of surgical techniques and continuously improve the quality of work in the operating theatre.

## New solutions


The design and production process is based on modern technologies that ensure safety and high quality of our products. Simultaneously, in order to provide products that outpace needs of our the most demanding clients, we are always in search for new solutions.

## Production process

The production process is maintained according to the highest standards, certified by the Quality Management System Certificates EN ISO 13485:2016. All our electrosurgical units are CE labelled.

We offer high quality and reliable electrosurgical units.

For confidence in performance.



# safety systems

## confidence through safety

To get through a successful operating procedure is always a matter of both - a medical team and instruments the surgeon can relay on. The rule is to pick out instruments in the same way as surgical team, by choosing those that ensure smooth running of any procedure and back-up in case of emergency.

AutoTest	EMED systems, after being switched on, perform a comprehensive internal test, whose result is displayed in the form of a precise, legible message.
Overload Protection	Supervising the temperature of the most critical components of the system makes it possible to avoid damaging the system, even if it is being heavily exploited.
Automatic operation Spectrum Result	Automatic regulation of the delivered power depending on the needs. A microprocessor monitors the tissue parameters on an ongoing basis to achieve the optimum cutting effect.
NEM System	The NEM system controls the quality of the contact of a return electrode with the patient's skin. Using the NEM system, as well as splitted neutral electrodes, has totally eliminated the problem of burning the patient's skin while applying a patient neutral electrode.
EndoDetect	Fully controlled resection process for polypectomy procedures. The system of snare fitting detection makes it impossible to activate the cutting mode when the snare does not fit tightly around the tissue. By preventing an accidental electrical current activation, Endo Detect minimizes the risk of perforation, ensuring maximum safety of the surgical procedures.
Defibrillator Proof	All EMED systems are resistant to the defibrillation impulse.
Constant surveillance of electrical parameters	A fast processor supervising the operation of the system provides an up-to-date assessment of the electrical parameters. Any possible disturbances are immediately reported.

# spectrum of use in electrosurgery

## result oriented performance



To facilitate the work with an electrosurgical system to the greatest extent possible, we have equipped spectrum with a number of functions to support the surgeon during the procedure.

- the **SmartDevice System** detects and identifies the connected instrument. It automatically adjusts the appropriate operating modes and output parameters to the connected instrument,
- the **SpectrumResult** solution maintains the output settings so as to obtain the desired effect regardless of the surgical conditions,
- ThermoStapler® II new generation of vessel sealing system,
- 10-inch **InTouch Screen** has integrated brightness adjustment and a choice of graphics versions of the screen,
- **universal SDS outputs** allow the use of monopolar and bipolar instruments in the same output,
- new specialist modes are adapted to specific applications in urology, arthroscopy and endoscopy,
- **simple software** updates via USB,
- communication with the user in different languages, **Voice Communication** commands inform about the operating status,
- **Spectrum Trolley** provides ergonomic work in the operating room.

**Electrosurgical unit SPECTRUM** is a state-of-the-art electrosurgical system that automatically adapts to the user's needs.

Owing to the solutions used in **SPECTRUM**, the user does not need to control power settings. The **SPECTRUM** system itself ensures that the output settings are maintained so as to obtain the desired result, regardless of the surgical conditions – **result oriented performance**.


Through active cooperation with our Customers we have created **SPECTRUM** – the first electrosurgical system which adapts to different surgical procedures so that each surgical intervention is as effective as possible.



# modes of operations

## setting performance goals

Get to know our electrosurgical operating modes to select them accordingly to your performance goals.



### MONO CUT


Monopolar cutting with different haemostasis effects. Effect 1 is used to cut the tissues when additional bleeding control is not needed. This cutting mode causes the least damage to tissues. Other effects include enhanced haemostasis levels. They are used when there is a need for more intensive bleeding control already in the course of cutting. The higher haemostasis levels control bleeding better, but exert a stronger thermal impact on the tissues. Instruments: monopolar electrodes, e.g. knife, loop or needle.

### MIXED CUT

Monopolar drying cutting. Alternating cutting and soft coagulation enable the cutting of severely bleeding tissues, while minimising blood loss. Instruments: monopolar electrodes, e.g. knife, loop or needle.

### PRECISE CUT

Precise monopolar cutting. Used when cutting small and precise structures. More gentle current enables higher cutting precision. Instruments: monopolar electrodes, e.g. knife, loop or needle.




### URO CUT

Monopolar cutting for urological procedures. The mode is used in the environment of non-conductive fluids, e.g. Purisol, glucose or distilled water. It is necessary for the TURP and TURB procedures. Instruments: monopolar urological resectoscope.

### HYSTERO CUT

Monopolar cutting for specialised gynaecological procedures (hysteroscopy). The mode is used in non-conductive fluids, e.g. Purisol or glucose. Instruments: monopolar hysteroscope - loop electrode.




### POLIPO CUT

Monopolar cutting for endoscopic procedures. Necessary for polyp removal. Alternating cutting and coagulation ensures optimum coagulation and reduces the risk of bleeding. Instruments: standard endoscopic loops, e.g. from Olympus or Pentax.

### PAPILLO CUT

Monopolar cutting for endoscopic procedures. Used for cutting Vater's papilla during a papillotomy procedure. It ensures safe cutting with the optimum hemostatic effect. Instrument: papillotome, e.g. from Olympus or Pentax.



### SOFT COAG

Low-voltage contact monopolar coagulation. This mode enables deep coagulation, reaching deeper than the other types. Instruments: monopolar electrodes, e.g. ball, spatula or lancet.

**SPRAY COAG**


High-voltage non-contact monopolar coagulation. It enables quick and effective coagulation of larger areas. It eliminates tissue adherence to the instrument. Instruments: monopolar electrodes, e.g. ball, spatula or lancet.

**HYBRID COAG**

Monopolar coagulation for contact and non-contact highvoltage applications. Instruments: monopolar electrodes, e.g. ball, spatula or lancet.

**FORCED COAG**

Contact monopolar coagulation. A traditional type of coagulation which enables quick and efficient coagulation of local bleeding. Instruments: monopolar electrodes, e.g. ball, spatula or lancet.

**ARGON COAG**


Argon-enhanced monopolar coagulation. The mode is used for non-contact coagulation of the surface of bleeding tissues. It eliminates smoke and smell. It ensures very shallow and gentle coagulation. Instruments: rigid argon electrodes for coagulation.

**PULSE ARGON**

Argon-enhanced pulsed monopolar coagulation. It is used in gastroenterology to control bleeding. It enables the precise delivery of energy doses exactly to the bleeding site. Instruments: flexible argon probes.

**ARGON CUT**


Argon-enhanced monopolar cutting. The use of argon reduces the amount of smoke and smell. The thermal damage to the tissues is reduced and bleeding control is improved. This function is particularly desirable during procedures that require intensive use of the unit. Instruments: needle or lancet type argon electrodes.

**BI-CUT**

Bipolar cutting with different effects of haemostasis. Special bipolar instruments are used for this mode. This mode is particularly recommended for procedures performed on neonates and patients with a heart pacemaker. Instruments: specialised standard instruments for bipolar cutting.

**URO BI-CUT**


Bipolar cutting for urological procedures TURP and TURB. Used for cutting and vaporisation of tissue. This mode is used in fluid environment. It requires the use of conductive fluids, e.g. normal saline or Ringer's solution. Instruments: bipolar urological resectoscope, loop or a vaporisation electrode.

**SOFT BI-COAG**

Low-voltage contact bipolar coagulation. In this mode, the current flows between the electrode tips and no passive electrode is required. Typically, it is used to seal single medium-sized blood vessels. Instruments: bipolar forceps, bipolar laparoscopic instruments.

# Argon

efficiency oriented performance




Argon in  
electrosurgery

Argon plasma coagulation ensures fast and efficient coagulation of large, heavily bleeding surfaces. Provides effectual devitalisation of tissues, e.g. neoplastic tumours. Argon plasma coagulation means less blood loss and less tissue damage. With penetration depth limited to 3 mm, it is particularly recommended for areas of high perforation risk.



Argon coagulation  
and cutting

EMED electro-surgical units are equipped with standard argon coagulation and pulse argon coagulation operating modes. Those modes allow to run any surgical intervention in safe and effective way. The option of adjustment of pulsation in argon coagulation enables more precision during endoscopic procedures. EMED electro-surgical units are also equipped with the option of monopolar cutting in argon shield.



Efficiency enhancing  
features

- immediate hemostasis helps efficiently coagulate large areas of bleeding surface
- penetration depth limited to approximately 3 mm minimizes risk of perforation
- tissue carbonization is minimal compared to standard electrocoagulation
- no tissue vaporization minimizes the risk of perforation
- no contact between the applicator and tissue means no tissue adhesion
- less surgical smoke gives good visibility of operating area
- reduced smoke eliminates unpleasant odors
- precise application of thermal energy results in reducing procedure time





Argon in  
endoscopy


The advantages of argon coagulation makes an ideal in endoscopic procedures, in gastroenterology in particular. Argon coagulation together with dedicated endoscopic modes of operation made by EMED allow safe and effective way of running surgical intervention.

For endoscopic interventions EMED dedicates particular flexible endoscopic probes. They are designed for multiple use and are available in a variety of sizes and lengths. Detailed list of argon instruments and accessories are available from current accessories catalogue.



Argon in  
open surgery

Argon coagulation is an excellent option for open surgery procedures, especially in case of glandular type of organs. For open surgery procedures we offer rigid argon applicators for coagulation and electrodes with needle or lancet ends for monopolar cutting in argon shield.



Electrosurgical units  
with argon

- Electrosurgical unit ATOM
- Electrosurgical unit ENDO
- Electrosurgical unit SPECTRUM
- Electrosurgical unit ES 350 with argon module
- Electrosurgical unit ES 350 with argon and ThermoStapler®


#### Applications

- gastroenterology
- general surgery
- pulmonology
- bronchoscopy
- otolaryngology

For more information about argon coagulation and its application see “Argon in electrosurgery” catalogue.

# ThermoStapler®

time oriented performance



ThermoStapler®

Due to the combination of the thermal and mechanic effects, ThermoStapler® offers a permanent alteration of collagen structure. Ensures higher resistance and closing. ThermoStapler® enables closing blood vessels with a diameter of 7 mm, however each time the user's assessment of the closing security is needed.



Time reducing features

- integrated system: argon coagulation and a vessel
- sealing module ensures cost reduction, there is no need to purchase any additional units
- high reliability of the system ensures maximum safety of work
- blood loss is significantly reduced
- the system allows to economize on sutures and staplers
- natural sealing of vessels
- multiple-use device, no limits on the number of uses
- no foreign body remains in the patient's body
- regulating the intensity of the effect enables the surgeon to select the most appropriate type of operation high reliability of the system ensures maximum safety of work



#### Advantages of ThermoStapler®

- physical pressure onto the vessel with the use of a bipolar clamps
- a generator deliver high frequency and low voltage current which
- does not cause sparks during the sealing process
- constant monitoring of the tissue impedance
- microprocessor control ensuring the repeated coagulation effect
- automatic switch-off of the generator at the moment of achieving
- the optimal sealing of the vessel, confirmed with a sound signal
- the possibility to select the level of intensity of the effect



#### ThermoStapler® instruments

EMED offers wide range of instruments dedicated to open surgery and laparoscopic procedures. All ThermoStapler® system instruments are multiple use and can be sterilized in an autoclave.

#### Electrosurgical units with ThermoStapler®

- Electrosurgical unit ATOM
- Electrosurgical unit SPECTRUM
- Electrosurgical unit ES 350 with argon and ThermoStapler®

#### Applications

- general surgery
- gynecology
- urology

For more information about vessel sealing system ThermoStapler® and its application see “ThermoStapler®” catalogue.

# Electrosurgical unit ATOM

compact but strong



SDS System  
 Touch screen  
 integrated argon  
 connection – SDSA  
 ThermoStapler® – vessel  
 sealing system  
 Automatic power regulation  
 Real-time measurements  
 of working parameters  
 Various configuration option  
 Bipolar resection  
 Special endoscopic modes  
 Power Monitor

Electrosurgical unit **ATOM** is a compact system which easily adapts to the user. It offers possibilities which have so far been provided by large and complicated electrosurgical systems only. Since it is possible to configure in any way the available operating modes, each user can set up an electrosurgical system that accurately meets the actual needs and requirements.

**ATOM** features monopolar and bipolar techniques in a variety of cutting and coagulation modes. The generator features highly-specialized operating modes, e.g. bipolar resection, an endoscopic cutting procedure and ThermoStapler® – a system for sealing large blood vessels. An integrated argon module brings the benefits from argon plasma coagulation and cutting in open, laparoscopic and endoscopic surgical procedures. **ATOM** is the first electrosurgical system which offers the advantages of both a compact size and a huge range of capabilities.



ATOM ensures an unprecedented ease of operation for the user with its:

- colour seven-inch display with a touch screen,
- innovative user interface,
- power monitor displaying average and real power,
- real-time monitoring of the output parameters of the unit,
- detection of the connected instruments (the SDS system),
- automatic selection of modes and settings to match the connected instrument,
- integrated argon connection,
- easy configuration of available modes and functions through an USB port,
- the possibility of an easy upgrade of its software to include new functions and operating modes.



# Electrosurgical unit ENDO

## for endoscopy



Touch screen  
 Advanced endoscopic modes  
 SpectrumResult  
 Monopolar modes  
 Bipolar mode  
 Argon coagulation  
 SDS System  
 NEM System  
 AUTOTEST

Electrosurgical unit ENDO is the only highly specialized for endoscopic procedures on the market. It has been attentively design to meet all requirements of endoscopic operation room. It enables to perform all endoscopic procedures that requires monopolar endoscopic cut, bipolar and monopolar coagulation. Integrated argon module allows to use all the advantages of argon coagulation in endoscopic procedures. SDS system recognizes

connected instruments and automatically selects the optimal working parameters. endo is equipped with touch screen display. Neutral electrode monitoring system [NEM] together with EMED SAFE disposable electrodes, guarantees maximum safety during procedures. endo is supported by specially design triple button foot switch for cut, coagulation and argon coagulation activation.



- advanced endoscopic cutting modes for polypectomy, sphincterectomy, mucosectomy
- specialized monopolar coagulation modes for endoscopic procedures
- continuous and pulsed endoscopic argon coagulation
- automatic regulation of all working parameters based on measurements in real time
- monitoring real and average power in endoscopic modes. It allows the operator to control unit parameters
- colour touch screen and SDS system for instruments detection
- neutral electrode monitoring system [NEM] and EMED SAFE electrodes guarantees maximum safety during procedure
- triple footswitch for independent activation of cutting, coagulation and argon plasma
- dedicated trolley with case for argon cylinder and handy basket for accessories and cables provides comfort and ergonomics in the operation room


#### ENDOSCOPIC CUTTING



#### ENDOSCOPIC COAGULATION



#### ARGON PLASMA



#### STANDARD MONOPOLAR MODES



# Electrosurgical unit ARTRO

## for arthroscopy



Touch screen  
 Advanced bipolar modes  
 for use in Saline solution  
 Monopolar modes  
 SDS System  
 NEM System  
 AUTOTEST  
 OVERLOAD  
 Automatic power regulation

Electrosurgical unit ARTRO is a highly specialized unit dedicated for bipolar and monopolar arthroscopic procedures. ARTRO allows vaporization, cutting and bipolar coagulation in Saline solution. It also offers monopolar cutting and coagulation modes which are necessary in some arthroscopic operations. Automatic regulation of output parameters guarantees fast and effective work. SDS system recognizes connected


instruments and automatically selects the optimal working parameters. The operating modes and the output settings are selected automatically to the connected instrument. ARTRO is equipped with touch screen display. Neutral electrode monitoring system [NEM] together with EMED SAFE disposable electrodes, guarantees maximum safety during procedures.





- advanced bipolar cutting and vaporization modes for arthroscopic procedures in Saline solution
- versatile due to monopolar cutting and coagulation modes
- automatic regulation of all working parameters based on measurements in real time
- colour touch screen provides quick and easy regulation of parameters
- instruments detection system facilitates the preparation for procedure
- after connecting SDS instrument the unit automatically selects suggested settings
- neutral electrode monitoring system [NEM] and EMED SAFE electrodes guarantees maximum safety during procedure selects suggested settings
- dedicated trolley equipped with shelves for accessories and materials necessary for the operation

BIPOLAR  
CUTTING



BIPOLAR  
VAPORIZATION



BIPOLAR  
COAGULATION




MONOPOLAR CUTTING  
AND COAGULATION



# Electrosurgical unit ATOM Smart

pure performance




Touch screen  
 SDS System  
 Compact size  
 Monopolar and bipolar  
 Excellent performance  
 Precise regulation  
 NEM system

Electrosurgical unit ATOM Smart is an electrosurgical device which uses high-frequency current flow for cutting and coagulating tissue. The unit allows to perform procedures in monopolar and bipolar modes. ATOM Smart electrosurgical system is controlled with a touch screen which is designed to provide the user with easy access to all functions. The settings or modes of operation are changed by touching icons on the screen. To ensure maximum

ease of operation there are no extra buttons or knobs. ATOM Smart is equipped with two universal SDS outputs, which identifies the connected instrument, automatically adjusts the appropriate operating modes and output parameters to the connected instrument. Neutral electrode monitoring system [NEM] together with EMED SAFE disposable electrodes, guarantees maximum safety during procedures.

# Electrosurgical unit ES 350

with argon and ThermoStapler®




ThermoStapler®  
 Argon  
 Bipolar resection  
 9 programs  
 Endoscopic cutting mode  
 MultiSwitch  
 SpectrumResult  
 NEM System  
 Overload Protection

Electrosurgical unit ES 350 with argon and ThermoStapler® is an electrosurgical system that meets the expectations of the most demanding medical professionals. It combines the benefits of electrosurgery, argon coagulation and a blood vessel closing system in one unit. It is equipped with two monopolar sockets, a bipolar socket and an independent ThermoStapler® socket. Enables operating with the use of various instruments without changing them during the procedure.

The unit follows an advanced operating modes such as bipolar resection. ES 350 customizes settings for different interventions by using 9 program storage positions. ES 350 with ARGON and Thermo-Stapler® system is equipped with SpectrumResult System for an automatic adjustment of power to maximizes the effectiveness of operation. Neutral Electrode Monitor System and Autotest ensures maximum safety during operation.

# Electrosurgical unit ES 350

with argon module



Argon

Bipolar resection

9 programs

Endoscopic cutting mode

MultiSwitch

SpectrumResult

NEM System

Overload Protection


An electrosurgical generator which meets all the requirements of a modern surgery.

The unit enables monopolar operation in standard modes and in ARGON enhanced cutting and coagulation. It is excellent both for open surgery and endoscopic and laparoscopic procedures. Equipped with SpectrumResult system for an automatic adjustment of power, it maximizes the

effectiveness of operation. The safety systems: Neutral Electrode Monitor and AutoTest guarantee security of performed procedures. The option of programming settings and modes additionally increases the comfort of work with **ES 350 with an ARGON module**. It enables remote program change via MultiSwitch.

# Electrosurgical unit ES 350

designed for modern surgery



Bipolar resection  
Endoscopic program  
MultiSwitch  
Spectrum Result  
NEM System  
Overload Protection

Electrosurgical unit ES 350 offers a variety of operating modes for open procedures, as well as endoscopy and laparoscopy. It enables monopolar and bipolar operation. It is also equipped with additional functions for specialized procedures in urology, arthroscopy and endoscopy. The unit is equipped with system that

enhances the effectiveness of operation: automatic power adjustment Spectrum Result. The neutral electrode monitoring NEM System and AutoTest guarantee maximum safety during surgical procedure. The generator enables remote change of programs via MultiSwitch.

# Electrosurgical unit ES 300

designed for outpatient clinics



Monopolar modes  
Bipolar mode  
Bipolar AutoStart  
Bipolar AutoStop  
Surgilogic  
PowerStart  
NEM System  
Overload Protection

A new version of ES 300 unit. It operates in monopolar cutting mode with different levels of haemostasis as well as mode of soft and forced coagulation. ES 300 enables bipolar coagulation controlled with footswitch or automatically with AutoStart or AutoStop function. It is equipped with PowerStart

system which facilitates the start of a procedure cutting and SURGILOGIC - automatic power adjustment. Supported by Neutral Electrode Monitor NEM, which controls the quality of neutral electrode application and AutoTest, an internal test conducted each time after power system is switched on.

# Electrosurgical unit ES 120

designed for small surgical procedures for outpatient clinics



Modern desing  
Automatic power regulation  
Excellent performance  
Monopolar modes  
Bipolar mode  
NEM System  
Overload Protection  
Footswitch recognition

The **ES 120** has separate control panel for operation in monopolar cut and coagulation, and bipolar coagulation mode. The unit is also equipped with a multiple footswitch. The generator operates in the following modes: pure cutting, cutting with haemostasis, monopolar soft and forced coagulation and bipolar coagulation. **ES 120**

is equipped with Power Start System, neutral electrode monitoring NEM, and AutoTest an internal test conducted each time after the power is switched on. Altogether with smart appearance and user-friendly applications makes **ES 120** an excellent choice for gynecological, dermatological and cosmetological surgeries, as well as for small surgical procedures for outpatient clinics.

# ARIA

## smoke evacuation system



Operating smoke is always an integral part of electrosurgical operation. Smoke generated during cutting and coagulation of tissues reduces the visibility of operating field. It also decreases the comfort of work in the operation room. Electrosurgical smoke, in a longer period of time, may have also negative influence on the health of people using electrosurgery in their daily work.

To improve comfort of operations and ensure the safety of the whole operating team it is recommended to use smoke evacuation system during electrosurgery. Aria system assures safe and comfortable work in the operation room. Aria is totally integrated with Spectrum electrosurgical system and it is also compatible the spectrum and other units EMED.



### Advantages


- quiet operation of the device guarantees no distractions during operation
- 35-hour filter life reduces the costs of procedures
- comprehensive three-port filter allows to connect 6mm, 9.5mm or 22mm accessories without using of additional adapters
- integration with Spectrum electrosurgical system allows automatic activation of suction while electrosurgical instrument is activated
- manual or via foot switch control available for other electrosurgical units
- three levels of air flow regulation
- easy to use





# WATERFALL

## endoscopic irrigation pump



Good visibility and clean operating site are of prime importance in all types of surgery.

In endoscopic procedures, the operation site may be contaminated with chime or blood residue in the gastrointestinal tract.

The possibility of irrigating the operating site is a very important element of any endoscopic procedure. The WATERFALL endoscopic irrigation pump was designed to ensure full comfort and safety during surgeries.

The pump makes it possible to quickly rinse the gastrointestinal tract using physiological saline or sterile water directly through the endoscope's rinsing channel or by using endoscopic instruments.

### Advantages


- good visibility and clean operating site
- quick and easy fluid flow adjustment from the pump control panel
- possibility of activation using the footswitch or from the control panel
- low-noise, does not cause nuisance in the procedure room
- small size and possibility of mounting of the device on the endo trolley, which gives access to the pump during the procedure

# Electrosurgical accessories


Select electrosurgical instruments and accessories to reach your performance goals. We offer assortment of accessories for electrosurgery, including:

monopolar accessories · bipolar accessories · laryngology · arthroscopy · ThermoStapler®  
vessel sealing clamps · neutral electrodes · monopolar and bipolar cables · argon coagulation instruments · papillotomy instruments · polypectomy instruments · other endoscopic and laparoscopic instruments

The full range of our accessories is presented in our **Accessories for Electrosurgery** catalogue. We also invite to purchase our products through emed official website: [www.emed.pl](http://www.emed.pl)



# contact us



## EMED SP. Z O. O. SP. K.

Ryżowa 69a,  
05-816 Opacz-Kolonia  
Poland  
tel: + 48 22 723 08 00  
export@emed.pl  
www.emed.pl



EMED products are available all over the world. See [www.emed.pl](http://www.emed.pl) for contact details.



Attention! This brochure does not replace instructions for use! Refer to instructions for use!

Copyright© EMED.

All rights reserved. Any copying, distribution, publishing in whole or in part without written EMED permission is prohibited.

REF NO 2021.11.K04.EN