



## Electrosurgical unit atom smart

Small size	<b>SAFETY FEATURES</b>	
Monopolar and bipolar modes	Protection against electric shock: Class/Grade	I/CF
Recognition of the footswitch	Low-frequency leakage currents	according to IEC 60601-1
Separate bipolar panel	High-frequency leakage currents	according to IEC 60601-2-2
Automatic power regulation	Generator operation frequency	333 [kHz]
	Defibrillation impulse resistance	according to EN 60601-1
	NEM System	neutral electrode control
	AUTOTEST	self-test of generator and accessories after switching the power on
	OVERLOAD	generator overload protection
SmartDevice System - detects and identifies the connected instrument	<b>CUTTING</b>	<b>MAX 200W</b>
	Pure monopolar cutting	9 effects
	Blend monopolar cutting	9 effects
Intuitive touch screen	<b>COAGULATION</b>	<b>MAX 200W</b>
Wireless footswitch	Soft monopolar coagulation	9 effects
Automatic power regulation	Forced monopolar coagulation	9 effects
Real time control of outputs parameters	<b>BIPOLAR COAGULATION</b>	<b>MAX 120W</b>
Easy software upgrade	Bipolar coagulation	9 effects
	<b>WEIGHT AND DIMENSIONS</b>	
	Length	385 [mm]
	Width	305 [mm]
	Height	141 [mm]
	Weight	5 [kg]
	<b>POWER SUPPLY</b>	
	Supply voltage	220-240 [V] +/- 10% or 110-120 [V] +/- 10%
	Rated power input	600 [VA]