

gastroenterology

safety and precision in endoscopic procedures



spectrum
result oriented performance

spectrum Spectrum is a state-of-the-art electrocautery system that automatically adapts to the user's needs so that each surgical intervention is as effective as possible.


SmartDevice System SmartDevice System detects and identifies the connected instrument. It automatically adjusts the appropriate operating modes and output parameters.

Endoscopic modes Modes available in EMED diathermies were specially designed to ensure high precision in endoscopic procedures due to fully controlled cutting process. Instant measurement of all cutting parameters guarantees efficient and safe operation in all conditions.

EMED Endo-Detect EndoDetect System minimises the risk of perforation by preventing accidental activation of current flow. It ensures the safety of the performed procedures.

ENDO ARGON This mode particularly recommended for areas of high perforation risk.


PULSE ARGON It is specially recommended for endoscopic surgery. Argon coagulation enables.



POLIPO CUT / PAPILOCUT

Indications:

- polypectomy
- sphincterotomy - biliary pancreatitis,
- endoscopic mucosal resection (EMR),
- endoscopic submucosal dissection (ESD),
- haemostasis of surface haemorrhage.



ENDO ARGON / PULSE ARGON

Indications:

- residual tissue removal and haemostasis of bleeding areas after polypectomy,
- bleeding from angiodysplastic lesions,
- stent ingrowth,
- bleeding from malignant lesions in the stomach and colon, colitis.

Endoscopic loop for polypectomy, oval, 2.3mm channel



E 17/25 big oval, Ø 40mm E 17/23 medium oval, Ø 30mm E 17/22 small oval, Ø 25mm

Endoscopic loop for polypectomy, hexagonal, 2.3mm channel



E 18/25 big hexagonal, Ø 40mm E 18/23 medium hexagonal, Ø 30mm E 18/22 small hexagonal, Ø 25mm



R1-100 Endoscopic loop handle, for E17-x, E18-x loops



E30-2 - Papillotomy snare loop 25mm, for R2-50 handle

E30-3 - Papillotomy snare loop 30mm



R2-50 Papillotomy loop handle




281-035 Monopolar endoscopic cable, L: 3m, 3mm female, SDS plug

Flexible argon probe, reusable



932-148 Ⓞ Ø 1.5 mm, length 1.5m 932-149 Ⓞ Ø 2.3 mm, length 2.2 m 932-150 Ⓞ Ø 3.2 mm, length 2.2 m 932-151 Ⓞ Ø 1.5 mm, length 3 m 932-152 Ⓞ Ø 2.3 mm, length 3 m



432-465 Monopolar cable for argon flexible electrode, L: 3.5m, SDS/LuerLock plug, flat connector

See our catalogue for a comprehensive range of accessories.
For more information, contact the regional representatives or visit www.emed.pl

EMED SP. Z O. O. SP. K. | ul. Ryżowa 69a | 05-816 Opacz Kolonia, Poland
tel: +48 (22) 723 08 00 | fax: +48 (22) 723 00 81 | email: export@emed.pl | www.emed.pl



100-013 Electrosurgical Unit spectrum



080-100 SpectrumLine trolley for EMED ESU, with case for argon cylinders (2x 5L/10L), basket



100-313 MultiSwitch, 2-button footswitch, wireless