


dermatology

a compact solution for specialist practice



ES120



Electrosurgical device

ES120 is an electrosurgical device which uses high-frequency current flow to cut and coagulate skin lesions.

Versatile use

The device is perfectly suited for dermatology, cosmetic procedures and dermatosurgical procedures.

Comfort of use

Small size of the device and its aesthetic appearance, combined with reliability and ease of use ensure the highest comfort of work.

Safety

The internal test of the device and connected accessories, as well as the neutral electrode connection control system guarantee maximum safety for the user and the patient.

Applications in dermatology:

- teleangiectasia
- unnecessary hair removal
- viral, seborrhoeic, flat verrucae
- molluscum contagiosum
- spider naevi
- fibromas
- atheromas
- miliums
- sample taking for tests
- other skin lesions

Advantages of the ES120 device:

- aesthetic design and small size
- effective and rapid operation
- precision
- low costs of use

accessories for ES120 device

dermatology equipment



541-800 Ball electrode, angled, 3mm, handle 1.6mm



540-600 Wire loop electrode, 8mm, handle 1.6mm



548-200 Spatula electrode straight, 2 x 6 mm,
L: 55mm, handle 1.6mm



541-600 Needle electrode, angled 30°, 1.7mm,
handle 1.6mm



328-030 Electrode handle 1.6mm, small,
no switches, cable 3m, 3-pin plug



605-021 Bipolar forceps, angled, 0.5mm, L: 110mm



non-stick



351-031 Bipolar cable, L: 3m, straight connector,
2-pin 29mm plug



100-301 One-pedal footswitch, for ES120, cable 5m,
6-pin plug



51301018 Handle type F for epilation needles 1.25mm,
4mm plug, separate cable 1.5m




F-type



- F3I Electrolysis needles, disposable, sterile, insulated, 1.25mm shank, 10 pcs., 0.003"
- F3 Electrolysis needles, disposable, sterile, non-insulated, 1.25mm shank, 10 pcs., 0.003"
- F3G Gold-plated telangiectasia needles; type F; 10 pcs.; 0.003"



294-030 Disposable neutral electrode cable, L: 3m,
EU 6.3mm plug



812-80H Disposable neutral electrode EMED SAFE,
hydrogel, split, for adults and children,
176x122mm, 110cm2, 10x5pcs.

See our catalogue for a comprehensive range of accessories.
For more information, contact the regional representatives or visit www.emed.pl